How to use bimaps from the ".db" annotation packages

Marc Carlson, Herve Pages, Seth Falcon, Nianhua Li

March 18, 2015

1 Introduction

1.0.1 Purpose

AnnotationDbi is used primarily to create mapping objects that allow easy access from R to underlying annotation databases. As such, it acts as the R interface for all the standard annotation packages. Underlying each AnnotationDbi supported annotation package is at least one (and often two) annotation databases. AnnotationDbi also provides schemas for theses databases. For each supported model organism, a standard gene centric database is maintained from public sources and is packaged up as an appropriate organism or "org" package.

1.0.2 Database Schemas

For developers, a lot of the benefits of having the information loaded into a real database will require some knowledge about the database schema. For this reason the schemas that were used in the creation of each database type are included in AnnotationDbi. The currently supported schemas are listed in the DBschemas directory of AnnotationDbi. But it is also possible to simply print out the schema that a package is currently using by using its "_dbschema" method.

There is one schema/database in each kind of package. These schemas specify which tables and indices will be present for each package of that type. The schema that a particular package is using is also listed when you type the name of the package as a function to obtain quality control information.

The code to make most kinds of the new database packages is also included in AnnotationDbi. Please see the vignette on SQLForge for more details on how to make additional database packages.

1.0.3 Internal schema Design of org packages

The current design of the organism packages is deliberately simple and gene centric. Each table in the database contains a unique kind of information and also an internal identifier called _id. The internal _id has no meaning outside of the context of a single database. But _id does connect all the data within a single database.

As an example if we wanted to connect the values in the genes table with the values in the kegg table, we could simply join the two tables using the internal _id column. It is very important to note however that _id does not have any absolute significance. That is, it has no meaning outside of the context of the database where it is used. It is tempting to think that an _id could have such significance because within a single database,

it looks and behaves similarly to an entrez gene ID. But _id is definitely NOT an entrez gene ID. The entrez gene IDs are in another table entirely, and can be connected to using the internal _id just like all the other meaningful information inside these databases. Each organism package is centered around one type of gene identifier. This identifier is found as the gene_id field in the genes table and is both the central ID for the database as well as the foreign key that chip packages should join to.

The chip packages are 'lightweight', and only contain information about the basic probe to gene mapping. You might wonder how such packages can provide access to all the other information that they do. This is possible because all the other data provided by chip packages comes from joins that are performed by AnnotationDbi behind the scenes at run time. All chip packages have a dependency on at least one organism package. The name of the organism package being depended on can be found by looking at its "ORGPKG" value. To learn about the schema from the appropriate organism package, you will need to look at the "_dbschema" method for that package. In the case of the chip packages, the gene_id that in these packages is mapped to the probe_ids, is used as a foreign key to the appropriate organism package.

Specialized packages like the packages for GO and KEGG, will have their own schemas but will also adhere to the use of an internal _id for joins between their tables. As with the organism packages, this _id is not suitable for use as a foreign key.

For a complete listing of the different schemas used by various packages, users can use the available.dbschemas function. This list will also tell you which model organisms are supported.

```
require(org.Hs.eg.db)
## Loading required package: org.Hs.eg.db
## Loading required package:
                              {\it Annotation Dbi}
## Loading required package: stats4
## Loading required package: BiocGenerics
## Loading required package: parallel
##
## Attaching package: 'BiocGenerics'
##
## The following objects are masked from 'package:parallel':
##
      clusterApply, clusterApplyLB, clusterCall, clusterEvalQ, clusterExport,
##
      clusterMap, parApply, parCapply, parLapply, parLapplyLB, parRapply,
##
##
     parSapply, parSapplyLB
##
## The following object is masked from 'package:stats':
##
##
      xtabs
##
## The following objects are masked from 'package:base':
##
      Filter, Find, Map, Position, Reduce, any Duplicated, append, as.data.frame,
##
##
      as.vector, cbind, colnames, do.call, duplicated, eval, evalq, get,
      intersect, is.unsorted, lapply, mapply, match, mget, order, paste, pmax,
##
      pmax.int, pmin, pmin.int, rank, rbind, rep.int, rownames, sapply, setdiff,
##
##
      sort, table, tapply, union, unique, unlist, unsplit
##
```

```
## Loading required package: Biobase
## Welcome to Bioconductor
##
      Vignettes contain introductory material; view with 'browseVignettes()'.
##
      cite Bioconductor, see 'citation("Biobase")', and for packages
##
      'citation("pkqname")'.
##
##
## Loading required package: GenomeInfoDb
## Loading required package: S4Vectors
## Loading required package: IRanges
##
## Attaching package: 'AnnotationDbi'
##
## The following object is masked from 'package:GenomeInfoDb':
##
##
      species
##
## Loading required package:
require(AnnotationForge)
## Loading required package: AnnotationForge
available.dbschemas()
```

2 Examples

2.0.4 Basic information

The AnnotationDbi package provides an interface to SQLite-based annotation packages. Each SQLite-based annotation package (identified by a ".db" suffix in the package name) contains a number of AnnDbBimap objects in place of the environment objects found in the old-style environment-based annotation packages. The API provided by AnnotationDbi allows you to treat the AnnDbBimap objects like environment instances. For example, the functions [[, get, mget, and ls all behave the same as they did with the older environment based annotation packages. In addition, new methods like [, toTable, subset and others provide some additional flexibility in accessing the annotation data.

```
library("hgu95av2.db")
##
```

The same basic set of objects is provided with the db packages:

```
## [9] "hgu95av2ENSEMBL"
                                 "hgu95av2ENSEMBL2PR0BE"
## [11] "hgu95av2ENTREZID"
                                 "hgu95av2ENZYME"
## [13] "hgu95av2ENZYME2PR0BE"
                                 "hgu95av2GENENAME"
## [15] "hgu95av2GO"
                                 "hgu95av2G02ALLPR0BES"
## [17] "hgu95av2G02PR0BE"
                                 "hgu95av2MAP"
## [19] "hgu95av2MAPCOUNTS"
                                 "hgu95av20MIM"
## [21] "hgu95av2ORGANISM"
                                 "hgu95av20RGPKG"
## [23] "hgu95av2PATH"
                                 "hgu95av2PATH2PR0BE"
## [25] "hgu95av2PFAM"
                                 "hgu95av2PMID"
## [27] "hgu95av2PMID2PR0BE"
                                 "hgu95av2PROSITE"
## [29] "hgu95av2REFSEQ"
                                 "hgu95av2SYMB0L"
## [31] "hgu95av2UNIGENE"
                                 "hgu95av2UNIPROT"
## [33] "hgu95av2_dbInfo"
                                 "hgu95av2_dbconn"
## [35] "hgu95av2_dbfile"
                                 "hgu95av2_dbschema"
```

Exercise 1

Start an R session and use the library function to load the hgu95av2.db software package. Use search() to see that an organism package was also loaded and then use the approriate "_dbschema" methods to the schema for the hgu95av2.db and org.Hs.eg.db packages.

It is possible to call the package name as a function to get some QC information about it.

```
qcdata = capture.output(hgu95av2())
head(qcdata, 20)
    [1] "Quality control information for hgu95av2:"
    [2] ""
    [3] ""
##
    [4] "This package has the following mappings:"
##
##
    [5] ""
    [6] "hgu95av2ACCNUM has 12625 mapped keys (of 12625 keys)"
##
    [7] "hgu95av2ALIAS2PROBE has 33755 mapped keys (of 103735 keys)"
    [8] "hgu95av2CHR has 11540 mapped keys (of 12625 keys)"
##
    [9] "hgu95av2CHRLENGTHS has 93 mapped keys (of 93 keys)"
## [10] "hgu95av2CHRLOC has 11474 mapped keys (of 12625 keys)"
## [11] "hgu95av2CHRLOCEND has 11474 mapped keys (of 12625 keys)"
## [12] "hgu95av2ENSEMBL has 11460 mapped keys (of 12625 keys)"
## [13] "hgu95av2ENSEMBL2PROBE has 9814 mapped keys (of 28553 keys)"
## [14] "hgu95av2ENTREZID has 11543 mapped keys (of 12625 keys)"
## [15] "hgu95av2ENZYME has 2125 mapped keys (of 12625 keys)"
## [16] "hgu95av2ENZYME2PROBE has 786 mapped keys (of 975 keys)"
## [17] "hgu95av2GENENAME has 11543 mapped keys (of 12625 keys)"
## [18] "hgu95av2GO has 11245 mapped keys (of 12625 keys)"
## [19] "hgu95av2GO2ALLPROBES has 17214 mapped keys (of 18826 keys)"
## [20] "hgu95av2G02PROBE has 12920 mapped keys (of 14714 keys)"
```

Alternatively, you can get similar information on how many items are in each of the provided maps by looking at the MAPCOUNTs:

hgu95av2MAPCOUNTS

To demonstrate the *environment* API, we'll start with a random sample of probe set IDs.

```
all_probes <- ls(hgu95av2ENTREZID)
length(all_probes)

## [1] 12625

set.seed(0xa1beef)
probes <- sample(all_probes, 5)
probes

## [1] "31882_at" "38780_at" "37033_s_at" "1702_at" "31610_at"</pre>
```

The usual ways of accessing annotation data are also available.

```
hgu95av2ENTREZID[[probes[1]]]

## [1] "9136"

hgu95av2ENTREZID$"31882_at"

## [1] "9136"

syms <- unlist(mget(probes, hgu95av2SYMBOL))

syms

## 31882_at 38780_at 37033_s_at 1702_at 31610_at

## "RRP9" "AKR1A1" "GPX1" "IL2RA" "PDZK1IP1"
```

The annotation packages provide a huge variety of information in each package. Some common types of information include gene symbols (SYMBOL), GO terms (GO), KEGG pathway IDs (KEGG), ENSEMBL IDs (ENSEMBL) and chromosome start and stop locations (CHRLOC and CHRLOCEND). Each mapping will have a manual page that you can read to describe the data in the mapping and where it came from.

?hgu95av2CHRLOC

Exercise 2

For the probes in 'probes' above, use the annotation mappings to find the chromosome start locations.

2.0.5 Manipulating Bimap Objects

Many filtering operations on the annotation *Bimap* objects require conversion of the *AnnDbBimap* into a *list*. In general, converting to lists will not be the most efficient way to filter the annotation data when using a SQLite-based package. Compare the following two examples for how you could get the 1st ten elements of the hgu95av2SYMBOL mapping. In the 1st case we have to get the entire mapping into list form, but in the second case we first subset the mapping object itself and this allows us to only convert the ten elements that we care about.

```
system.time(as.list(hgu95av2SYMBOL)[1:10])
## vs:
```

```
system.time(as.list(hgu95av2SYMBOL[1:10]))
```

There are many different kinds of *Bimap* objects in AnnotationDbi, but most of them are of class *AnnDbBimap*. All /RclassBimap objects represent data as a set of left and right keys. The typical usage of these mappings is to search for right keys that match a set of left keys that have been supplied by the user. But sometimes it is also convenient to go in the opposite direction.

The annotation packages provide many reverse maps as objects in the package name space for backwards compatibility, but the reverse mappings of almost any map is also available using revmap. Since the data are stored as tables, no extra disk space is needed to provide reverse mappings.

So now that you know about the revmap function you might try something like this:

```
as.list(revmap(hgu95av2PATH)["00300"])
## $`00300`
## [1] "35870_at" "36132_at"
```

Note that in the case of the PATH map, we don't need to use revmap(x) because hgu95av2.db already provides the PATH2PROBE map:

```
x <- hgu95av2PATH
## except for the name, this is exactly revmap(x)
revx <- hgu95av2PATH2PROBE
revx2 <- revmap(x, objName="PATH2PROBE")
revx2
## PATH2PROBE map for chip hgu95av2 (object of class "ProbeAnnDbBimap")
identical(revx, revx2)
## [1] TRUE
as.list(revx["00300"])
## $`00300`
## [1] "35870_at" "36132_at"</pre>
```

Note that most maps are reversible with revmap, but some (such as the more complex GO mappings), are not. Why is this? Because to reverse a mapping means that there has to be a "value" that will always become the "key" on the newly reversed map. And GO mappings have several distinct possibilities to choose from (GO ID, Evidence code or Ontology). In non-reversible cases like this, AnnotationDbi will usually provide a pre-defined reverse map. That way, you will always know what you are getting when you call revmap

While we are on the subject of GO and GO mappings, there are a series of special methods for GO mappings that can be called to find out details about these IDs. Term,GOID, Ontology, Definition,Synonym, and Secondary are all useful ways of getting additional information about a particular GO ID. For example:

```
Term("GD:0000018")
## Loading required package: GD.db
```

Exercise 3

Given the following set of RefSeq IDs: c("NG_005114","NG_007432","NG_008063"), Find the Entrez Gene IDs that would correspond to those. Then find the GO terms that are associated with those entrez gene IDs. org.Hs.eg.db packages.

2.0.6 The Contents and Structure of Bimap Objects

Sometimes you may want to display or subset elements from an individual map. A *Bimap* interface is available to access the data in table (*data.frame*) format using [and toTable.

```
head(toTable(hgu95av2G0[probes]))
##
     probe_id
                   go_id Evidence Ontology
## 1 1702_at GO:0006915
                              TAS
## 2 1702_at GO:0006955
                                         BP
                              TAS
## 3 1702_at GO:0007166
                              TAS
                                         BP
## 4 1702_at GD:0008283
                              TAS
                                         BP
## 5 1702_at GO:0038110
                              TAS
                                         BP
## 6 1702_at GD:0002437
                              TF.A
                                         BP
```

The toTable function will display all of the information in a *Bimap*. This includes both the left and right values along with any other attributes that might be attached to those values. The left and right keys of the *Bimap* can be extracted using Lkeys and Rkeys. If is is necessary to only display information that is directly associated with the left to right links in a *Bimap*, then the links function can be used. The links returns a data frame with one row for each link in the bimap that it is applied to. It only reports the left and right keys along with any attributes that are attached to the edge between these two values.

Note that the order of the cols returned by toTable does not depend on the direction of the map. We refer to it as an 'undirected method':

```
toTable(x)[1:6, ]
    probe_id path_id
## 1 1000_at
                04010
## 2 1000_at
                04012
## 3 1000_at
                04062
## 4 1000_at
                04114
## 5 1000_at
                04150
## 6 1000_at
                04270
toTable(revx)[1:6, ]
    probe_id path_id
```

```
## 1 1000_at 04010
## 2 1000_at 04012
## 3 1000_at 04062
## 4 1000_at 04114
## 5 1000_at 04150
## 6 1000_at 04270
```

Notice however that the Lkeys are always on the left (1st col), the Rkeys always in the 2nd col

For length() and keys(), the result does depend on the direction, hence we refer to these as 'directed methods':

```
length(x)
## [1] 12625
length(revx)
## [1] 229
allProbeSetIds <- keys(x)
allKEGGIds <- keys(revx)</pre>
```

There are more 'undirected' methods listed below:

Notice how they give the same result for x and revmap(x)

You might be tempted to think that Lkeys and Llength will tell you all that you want to know about the left keys. But things are more complex than this, because not all keys are mapped. Often, you will only want to know about the keys that are mapped (ie. the ones that have a corresponding Rkey). To learn this you want to use the mappedkeys or the undirected variants mappedLkeys and mappedRkeys. Similarly, the count.mappedkeys, count.mappedLkeys and count.mappedRkeys methods are very fast ways to determine how many keys are mapped. Accessing keys like this is usually very fast and so it can be a decent strategy to subset the mapping by 1st using the mapped keys that you want to find.

```
x = hgu95av2ENTREZID[1:10]
## Directed methods
mappedkeys(x)  # mapped keys

## [1] "1000_at" "1001_at" "1002_f_at" "1003_s_at" "1004_at"

## [6] "1005_at" "1006_at" "1008_f_at" "1009_at"

count.mappedkeys(x)  # nb of mapped keys

## [1] 9
```

```
## Undirected methods
mappedLkeys(x)  # mapped left keys

## [1] "1000_at" "1001_at" "1002_f_at" "1003_s_at" "1004_at"

## [6] "1005_at" "1006_at" "1008_f_at" "1009_at"

count.mappedLkeys(x)  # nb of mapped Lkeys

## [1] 9
```

If you want to find keys that are not mapped to anything, you might want to use isNA.

```
y = hgu95av2ENTREZID[isNA(hgu95av2ENTREZID)]  # usage like is.na()
Lkeys(y)[1:4]
## [1] "1007_s_at" "1047_s_at" "1089_i_at" "108_g_at"
```

Exercise 4

How many probesets do not have a GO mapping for the hgu95av2.db package? How many have no mapping? Find a probeset that has a GO mapping. Now look at the GO mappings for this probeset in table form.

2.0.7 Some specific examples

Lets use what we have learned to get information about the probes that are are not assigned to a chromosome:

```
x <- hgu95av2CHR
Rkeys(x)
## [1] "19" "12" "8" "14" "3" "2" "17" "16" "9" "X" "6" "1" "7"
## [14] "10" "11" "22" "5" "18" "15" "Y" "20" "21" "4" "13" "MT" "Un"
chroms <- Rkeys(x)[23:24]
chroms
## [1] "4" "13"
Rkeys(x) <- chroms</pre>
toTable(x)
##
        probe_id chromosome
## 1
       1029_s_at
## 2
       1036_at
                          4
## 3
         1058_at
                         13
## 4
         1065_at
                         13
## 5
         1115_at
                          4
## 6
         1189_at
                         13
## 7
         1198_at
                         13
## 8
        1219_at
                          4
## 9
      1220_g_at
                          4
## 10
        1249_at
                          4
                          4
## 11
         1285_at
## 12
        1303_at
```

```
## 13
      1325_at
## 14
       1348_s_at
                          13
## 15
       1369_s_at
                           4
## 16
        1377_at
                           4
## 17
       1378_g_at
                           4
## 18
       1451_s_at
                          13
## 19
        1503_at
                          13
## 20
       1507_s_at
                          4
## 21
       1527_s_at
                          13
## 22
        1528_at
                          13
## 23
         1529_at
                          13
## 24
       1530_g_at
                          13
## 25
        1531_at
                          13
## 26
       1532_g_at
                          13
## 27
       1538_s_at
                          4
## 28
        1542_at
                          4
## 29
       1545_g_at
                          13
## 30
        1567_at
                          13
## 31
       1570_f_at
                          13
## 32
       1571_f_at
                          13
## 33
        1593_at
                          4
## 34
         1597_at
                          13
## 35
       1598_g_at
                          13
## 36
         159_at
                           4
## 37
        1600_at
                           4
## 38
         1604_at
                           4
## 39
                           4
       1605_g_at
         1616_at
## 40
                          13
## 41
         1624_at
                           4
       1629_s_at
                          4
## 42
## 43
        1670_at
                          13
## 44
       1672_f_at
                          13
## 45
        1679_at
                          4
## 46
         1708_at
                           4
## 47
       1709_g_at
                           4
## 48
         170_at
                          13
## 49
         1720_at
                           4
## 50
       1721_g_at
## 51
       1731_at
                           4
## 52
         1732_at
                           4
## 53
         1819_at
                          13
## 54
       1828_s_at
                           4
## 55
        1836_at
                           4
## 56
       1883_s_at
                           4
## 57
       1888_s_at
                           4
## 58
        1900_at
                          13
## 59
       1905_s_at
                          13
```

## 60				
## 62				
## 63				
## 64				
## 65			1934_s_at	4
## 66			1943_at	4
## 67	##	65	1954_at	4
## 68	##	66	1963_at	13
## 69	##	67	1964_g_at	13
## 70	##	68	1987_at	4
## 71 1990_g_at	##	69	1988_at	4
## 72	##	70	1989_at	13
## 73	##	71	1990_g_at	13
## 74 2092_s_at	##	72	2044_s_at	13
## 75	##	73	2062_at	4
## 75	##	74		4
## 76				4
## 77				
## 78			_	
## 79				
## 80			_	
## 81				
## 82 31320_at 13 ## 83 31333_at 4 ## 84 31345_at 4 ## 85 31349_at 4 ## 86 31356_at 4 ## 87 31382_f_at 4 ## 89 31404_at 13 ## 90 31464_at 13 ## 91 31465_g_at 13 ## 92 31516_f_at 13 ## 93 31543_at 4 ## 94 31562_at 13 ## 95 31584_at 13 ## 96 31628_at 13 ## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 83				
## 84				
## 85				
## 86				
## 87 31382_f_at 4 ## 88 31404_at 13 ## 89 31408_at 4 ## 90 31464_at 13 ## 91 31465_g_at 13 ## 92 31516_f_at 13 ## 93 31543_at 4 ## 94 31562_at 13 ## 95 31584_at 13 ## 96 31628_at 13 ## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 88				
## 89				
## 90				
## 91 31465_g_at 13 ## 92 31516_f_at 13 ## 93 31543_at 4 ## 94 31562_at 13 ## 95 31584_at 13 ## 96 31628_at 13 ## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 92 31516_f_at 13 ## 93 31543_at 4 ## 94 31562_at 13 ## 95 31584_at 13 ## 96 31628_at 13 ## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 93			_	
## 94 31562_at 13 ## 95 31584_at 13 ## 96 31628_at 13 ## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 95				
## 96 31628_at 13 ## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 97 31631_f_at 4 ## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 98 31639_f_at 13 ## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 99 31640_r_at 13 ## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 100 31670_s_at 4 ## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 101 31684_at 4 ## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 102 31686_at 4 ## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 103 31706_at 4 ## 104 31744_at 4 ## 105 31753_at 13				
## 104 31744_at 4 ## 105 31753_at 13				
## 105 31753_at 13				
## 100 31/90_at 13				
	##	100	31/90_at	13

```
## 107
      31792_at
        31805_at
## 108
                           4
## 109 31811_r_at
                           4
## 110
       31847_at
                          13
## 111
       31849_at
                          13
## 112
       31851_at
                          13
## 113 31876_r_at
                           4
## 114
         31894_at
                           4
## 115 31969_i_at
                           4
## 116 31970_r_at
                           4
## 117 32006_r_at
                           4
## 118 32026_s_at
                           4
                           4
## 119
        32080_at
## 120
        32102_at
                          13
## 121
      32145_at
                          4
## 122 32146_s_at
                          4
## 123
        32147_at
                          13
## 124
        32148_at
                          13
## 125 32163_f_at
                          4
                           4
## 126 32180_s_at
## 127
        32220_at
                          13
## 128
       32299_at
                           4
                           4
## 129
       32349_at
## 130
      32353_at
                           4
## 131
                          4
        32357_at
## 132
        32368_at
                          13
## 133 32393_s_at
                          4
        32439_at
## 134
                          13
## 135
        32446_at
                          4
## 136
        32449_at
                           4
## 137
        32465_at
                           4
## 138
        32482_at
                          13
## 139
        32506_at
                           4
## 140
        32507_at
                           4
## 141
        32570_at
                           4
                           4
## 142
        32580_at
## 143
                           4
        32595_at
## 144
        32602_at
                           4
## 145
        32641_at
                          13
## 146
        32675_at
                          4
                           4
## 147
        32703_at
## 148
        32768_at
                          13
                           4
## 149
        32769_at
        32770_at
                           4
## 150
## 151
        32771_at
                           4
## 152
        32812_at
                           4
                           4
## 153
        32822_at
```

```
## 154
         32832_at
                           4
        32862_at
## 155
                           13
## 156
         32906_at
                           13
## 157
         32979_at
                           4
## 158 32986_s_at
                           13
## 159
         32998_at
                            4
## 160
         33013_at
                            4
                            4
## 161
         33050_at
## 162 33068_f_at
                            4
## 163 33069_f_at
                            4
## 164
         33100_at
                            4
         33150_at
                            4
## 165
## 166 33151_s_at
                            4
## 167
         33155_at
                            4
## 168
         33156_at
                            4
## 169
         33168_at
                           13
## 170 33171_s_at
                            4
## 171
         33172_at
                            4
## 172 33173_g_at
                            4
## 173
         33199_at
                           13
## 174
         33208_at
                           13
## 175
        33241_at
                            4
                            4
## 176
         33249_at
## 177
        33267_at
                            4
## 178
         33276_at
                           13
## 179
        33299_at
                           4
## 180
                           13
         33318_at
                            4
## 181
         33356_at
## 182
         33359_at
                            4
                            4
## 183
         33369_at
## 184 33370_r_at
                            4
## 185
         33382_at
                            4
## 186
         33483_at
                            4
## 187
         33488_at
                            4
## 188
         33490_at
                            4
                            4
## 189
         33494_at
## 190
                            4
         33519_at
## 191
         33520_at
                           13
## 192
         33525_at
                            4
## 193
                            4
         33526_at
                            4
## 194
         33529_at
## 195
         33536_at
                            4
## 196
                            4
         33544_at
## 197
         33564_at
                            4
## 198
         33576_at
                           13
## 199
         33584_at
                            4
## 200
                            4
         33596_at
```

```
## 201
        33657_at
## 202 33672_f_at
                          4
## 203 33673_r_at
                          4
## 204
        33687_at
                         13
## 205
        33700_at
                         13
## 206
       33733_at
                          4
## 207
        33791_at
                         13
## 208
       33823_at
                          4
## 209
       33827_at
                         13
## 210
       33837_at
                          4
## 211
      33859_at
                         13
## 212
                          4
       33975_at
## 213
        33990_at
                          4
## 214 33991_g_at
                          4
## 215
        33992_at
## 216
        33997_at
                          4
## 217
       34021_at
                          4
## 218
       34022_at
                          4
## 219
        34026_at
                         13
                          4
## 220
        34029_at
## 221
        34048_at
                          4
## 222
       34051_at
                         13
## 223
        34058_at
                          4
## 224
       34075_at
                          4
## 225
                          4
       34122_at
## 226
       34131_at
                          4
## 227
       34144_at
                          4
## 228
        34145_at
                          4
## 229
       34149_at
                          4
## 230 34170_s_at
                          4
## 231
        34181_at
                          4
## 232
        34198_at
                          4
## 233
       34211_at
                         13
## 234
        34239_at
                         13
## 235 34240_s_at
                         13
## 236
                          4
        34247_at
## 237
                          4
        34248_at
## 238 34275_s_at
                          4
## 239
        34284_at
                         13
## 240
       34307_at
                         13
                          4
## 241
      34319_at
## 242 34324_at
                         13
## 243
      34334_at
                         13
## 244
        34335_at
                         13
## 245
        34341_at
                          4
## 246 34342_s_at
                          4
## 247
                          4
        34353_at
```

```
## 248
         34398_at
                           13
## 249
         34411_at
                            4
## 250
         34423_at
                            4
## 251
         34459_at
                           13
## 252 34476_r_at
                            4
## 253
         34482_at
                            4
## 254
         34512_at
                            4
## 255
        34551_at
                            4
## 256
        34564_at
                            4
## 257
                            4
        34565_at
## 258
       34578_at
                           13
## 259
                           13
        34583_at
## 260
                            4
         34596_at
## 261 34637_f_at
                            4
## 262 34638_r_at
                           4
## 263
         34657_at
                           13
## 264
         34672_at
                           13
## 265
         34745_at
                           4
## 266
        34803_at
                           13
                            4
## 267
         34898_at
                            4
## 268 34953_i_at
## 269 34954_r_at
                            4
## 270
         34955_at
                           13
## 271
         34973_at
                            4
## 272
                            4
        34984_at
## 273
        34988_at
                            4
## 274
                            4
         35020_at
## 275
         35021_at
                            4
## 276
         35025_at
                            4
                            4
## 277
         35028_at
## 278
         35039_at
                            4
## 279
         35053_at
                            4
## 280
         35061_at
                            4
## 281
         35063_at
                            4
## 282
         35081_at
                           13
## 283
         35105_at
                           13
## 284
         35107_at
                           13
## 285
         35110_at
                           13
## 286
         35131_at
                           4
## 287
                           4
         35134_at
## 288
         35140_at
                           13
## 289
         35147_at
                           13
## 290
                            4
         35164_at
## 291
         35181_at
                            4
## 292 35182_f_at
                            4
## 293
         35193_at
                           13
## 294
         35213_at
                           13
```

```
## 295
         35214_at
## 296
         35215_at
                            4
## 297
         35220_at
                            4
## 298
         35285_at
                            4
## 299
         35306_at
                            4
## 300
         35344_at
                           13
## 301
         35356_at
                            4
## 302
         35357_at
                            4
## 303
         35371_at
                            4
## 304 35372_r_at
                            4
## 305
         35400_at
                           13
## 306
         35410_at
                            4
## 307 35435_s_at
                            4
## 308
         35437_at
                           4
## 309
         35469_at
                           13
## 310
         35470_at
                           13
## 311 35471_g_at
                           13
## 312
         35481_at
                           13
## 313
         35507_at
                           4
                            4
## 314
         35523_at
## 315 35554_f_at
                           13
## 316 35555_r_at
                           13
## 317
         35591_at
                            4
## 318
         35656_at
                           13
## 319
                            4
        35662_at
## 320
        35664_at
                            4
## 321
                            4
         35678_at
## 322
         35689_at
                            4
## 323
        35698_at
                            4
## 324
        35725_at
                           13
## 325
        35730_at
                            4
## 326
         35777_at
                            4
## 327
         35793_at
                            4
## 328
         35827_at
                            4
## 329
         35837_at
                            4
## 330
                            4
         35845_at
## 331 35871_s_at
                            4
## 332
         35877_at
                           13
## 333
         35904_at
                           13
## 334 35939_s_at
                           13
## 335
         35940_at
                           13
## 336
        35949_at
                           13
## 337
                           13
         35972_at
## 338
        35989_at
                            4
## 339
         35991_at
                            4
## 340
         36012_at
                           13
## 341
                            4
         36013_at
```

```
## 342
         36017_at
                           13
## 343
         36021_at
                            4
## 344
         36031_at
                           13
## 345
         36046_at
                            4
## 346
         36047_at
                            4
## 347
         36065_at
                            4
## 348
         36080_at
                            4
                            4
## 349
         36143_at
## 350
         36157_at
                            4
## 351
         36188_at
                           13
## 352
         36194_at
                            4
## 353
                           13
         36212_at
## 354
                            4
         36243_at
## 355 36247_f_at
                            4
## 356
         36269_at
                            4
## 357
         36274_at
                           13
## 358
         36358_at
                            4
## 359
         36363_at
                            4
## 360
         36433_at
                            4
                            4
## 361 36434_r_at
## 362
         36510_at
                           13
## 363
         36521_at
                           13
                            4
## 364
         36606_at
## 365
         36622_at
                            4
## 366
                            4
         36627_at
## 367
         36659_at
                           13
## 368
         36717_at
                            4
## 369
         36788_at
                           13
## 370
           367_at
                           13
                            4
## 371
         36814_at
## 372
         36830_at
                           13
## 373
         36913_at
                            4
## 374
         36914_at
                            4
## 375
         36915_at
                            4
                            4
## 376
         36918_at
## 377
                            4
         36939_at
                           13
## 378 36968_s_at
## 379
                            4
         36990_at
## 380
         37006_at
                            4
## 381
                            4
         37019_at
## 382
         37023_at
                           13
## 383
         37056_at
                            4
## 384
                            4
         37058_at
## 385
         37062_at
                            4
## 386
         37067_at
                           13
## 387
         37079_at
                           13
## 388
         37099_at
                           13
```

```
## 389
         37109_at
                          13
## 390
         37154_at
                           13
## 391
         37170_at
                           4
## 392
         37172_at
                           13
## 393
         37173_at
                           4
## 394
         37187_at
                            4
## 395
         37206_at
                           4
## 396
        37219_at
                           4
## 397
         37223_at
                           4
## 398
        37243_at
                           4
## 399
        37244_at
                           13
## 400
                           4
        37280_at
## 401
         37282_at
                           4
## 402 37291_r_at
                           4
## 403
         37303_at
                          13
## 404 37322_s_at
                           4
## 405 37323_r_at
                           4
## 406 37356_r_at
                           4
## 407
         37366_at
                           4
                           4
## 408
        37404_at
## 409
        37416_at
                           4
## 410
       37472_at
                           4
## 411
         37518_at
                          13
       37520_at
## 412
                           4
## 413 37521_s_at
                           4
## 414 37522_r_at
                           4
## 415
                          13
        37571_at
## 416
         37578_at
                           4
## 417
         37593_at
                          13
## 418
         37619_at
                           4
## 419
         37658_at
                          13
## 420 37707_i_at
                           4
## 421 37708_r_at
                           4
## 422
         37723_at
                           4
## 423
         37747_at
                           4
       37748_at
                           4
## 424
## 425
                           4
         37752_at
## 426
         37757_at
                           13
## 427
         37767_at
                           4
## 428
         37840_at
                           4
                           4
## 429
         37852_at
## 430
        37926_at
                           13
## 431
         37930_at
                           13
## 432
         37964_at
                           4
## 433
         38008_at
                           4
## 434
         38016_at
                            4
## 435
                            4
         38024_at
```

```
## 436 38025_r_at
                          4
## 437
        38035_at
                          13
## 438
        38065_at
                           4
## 439 38102_at
                          13
## 440
       38120_at
                           4
## 441
        38168_at
                           4
## 442
        38254_at
                           4
## 443 38304_r_at
                          13
## 444
       38353_at
                          13
## 445
        38375_at
                          13
## 446
        38438_at
                           4
## 447
                           4
        38485_at
                           4
## 448 38488_s_at
## 449
        38489_at
                           4
## 450
        38587_at
                           4
## 451
        38606_at
                           4
## 452
       38615_at
                          13
## 453
        38643_at
                           4
## 454
        38649_at
                          13
                           4
## 455
        38714_at
                           4
## 456
         38715_at
## 457
         38736_at
                           4
                           4
## 458 38751_i_at
## 459 38752_r_at
                           4
## 460
                           4
        38767_at
## 461
        38768_at
                           4
## 462
                           4
        38778_at
## 463
                           4
        38821_at
## 464
        38825_at
                           4
## 465
        38838_at
                           4
## 466
        38854_at
                           4
## 467
        38891_at
                           4
## 468
        38957_at
                          13
## 469
        38972_at
                          13
## 470
        38988_at
                           4
## 471
        39028_at
                          13
## 472
                          13
        39032_at
## 473
        39037_at
                           4
## 474
        39056_at
                           4
## 475
                           4
        39083_at
## 476
        39131_at
                          13
## 477
         39132_at
                           4
                           4
## 478 39208_i_at
## 479 39209_r_at
                           4
## 480
         39224_at
                           4
## 481
        39256_at
                          13
## 482
        39257_at
                          13
```

##	483	39269_at	13
##	484	39295_s_at	4
##	485	39297_at	13
##	486	39333_at	13
##	487	39337_at	4
##	488	39355_at	4
##	489	39369_at	4
	490	39380_at	4
	491	39382_at	4
	492	39405_at	13
		39469_s_at	13
	494	39475_at	4
	495		4
		39481_at	
	496	39488_at	13
		39489_g_at	13
	498	39535_at	4
	499	39536_at	4
	500	39554_at	4
	501	39555_at	4
##	502	39576_at	4
##	503	39579_at	13
##	504	39600_at	4
##	505	39634_at	4
##	506	39662_s_at	4
##	507	39665_at	4
##	508	39680_at	4
	509	39690_at	4
	510	39698_at	4
	511	39734_at	4
	512	39746_at	4
	513	39748_at	13
		39758_f_at	13
	515		
	516	39786_at	4
	517	39847_at	4
	518	39850_at	4
	519	39851_at	4
	520	39852_at	13
	521	39878_at	13
	522	39897_at	4
##	523	39924_at	13
##	524	39929_at	4
##	525	39955_at	13
##	526	39960_at	4
	527	39979_at	13
	528	40018_at	13
		40058_s_at	4
	_		

```
## 530 40059_r_at
                          4
## 531 40060_r_at
                          4
## 532
        40067_at
                          13
## 533
       40072_at
                          13
## 534
       40082_at
                          4
## 535
        400_at
                          13
## 536
        40114_at
                          4
## 537
       40121_at
                          4
## 538
        40148_at
                          4
## 539
        40180_at
                         13
## 540 40181_f_at
                         13
                          4
## 541
        40199_at
                          4
## 542 40217_s_at
## 543
        40218_at
                          4
## 544
       40225_at
## 545
       40226_at
                          4
## 546
      40272_at
                          4
## 547
      40310_at
                          4
## 548
       40312_at
                         13
## 549
       40323_at
                          4
## 550
       40349_at
                          4
## 551
        40354_at
                          4
## 552
        40392_at
                         13
## 553 40404_s_at
                         13
## 554
        40449_at
                          4
## 555
        40454_at
                          4
## 556
                          4
        40456_at
## 557
        40473_at
                          13
## 558
                          4
       40492_at
## 559
        40530_at
                          4
## 560
        40570_at
                         13
## 561 40576_f_at
                          4
## 562
        40633_at
                         13
## 563
        40681_at
                         13
                          4
## 564
        40697_at
                          4
## 565
       40710_at
                           4
## 566
        40711_at
## 567
        40727_at
                           4
## 568
        40746_at
                           4
                           4
## 569 40770_f_at
                           4
## 570
       40772_at
## 571
      40773_at
                          4
                          4
## 572
      40818_at
## 573
       40828_at
                          13
## 574
        40839_at
                          13
## 575
        40853_at
                          4
                           4
## 576 40880_r_at
```

```
## 577
        40893_at
                         13
## 578
          408_at
                         4
## 579 40908_r_at
                         13
## 580
       40943_at
                         4
## 581
      40970_at
                         13
## 582 40989_at
                          4
## 583
       40990_at
                          4
## 584
      40991_at
## 585 40992_s_at
                          4
## 586 40993_r_at
                          4
## 587 41014_s_at
                          4
                          4
## 588 41024_f_at
## 589 41025_r_at
                          4
## 590 41026_f_at
                         4
## 591
      41069_at
                         13
## 592
       41071_at
                         4
## 593 41104_at
                         4
## 594
      41118_at
                         13
## 595 41119_f_at
                         13
                         4
## 596 41145_at
                         4
## 597
      41148_at
## 598 41182_at
                         13
## 599 41191_at
                         4
## 600
      41276_at
                         13
## 601
      41277_at
                         13
## 602 41300_s_at
                         13
## 603
       41301_at
                         13
                         4
## 604
       41308_at
## 605 41309_g_at
                         4
## 606
        41317_at
                         13
## 607 41318_g_at
                         13
## 608
        41319_at
                         13
## 609 41376_i_at
                         4
## 610 41377_f_at
                          4
## 611
       41391_at
                          4
                          4
## 612 41392_at
## 613 41402_at
                          4
## 614 41434_at
                          4
## 615 41436_at
                         13
## 616 41456_at
                         4
## 617
      41459_at
                         13
## 618
       41470_at
                         4
## 619 41491_s_at
                         13
## 620 41492_r_at
                         13
## 621 41493_at
                         13
## 622
      41534_at
                          4
## 623 41555_at
                          4
```

```
## 624 41556_s_at
                            4
## 625
                            4
         41585_at
## 626 41667_s_at
                            13
## 627 41668_r_at
                            13
## 628
         41697_at
                            4
## 629
         41801_at
                            4
## 630
         41806_at
                            4
## 631
         41860_at
                           13
## 632
           431_at
                            4
## 633
           504_at
                            4
## 634
         507_s_at
                            4
                             4
## 635
           579_at
## 636
           618_at
                             4
## 637
           630_at
                             4
## 638
         631_g_at
                             4
## 639
           655_at
                             4
## 640
         690_s_at
                             4
## 641
         692_s_at
                             4
## 642
         764_s_at
                             4
## 643
           820_at
                            4
## 644
           886_at
                            4
## 645
           931_at
                           13
## 646
         936_s_at
                            4
## 647
         948_s_at
                            4
## 648
                            13
           963_at
## 649
           975_at
                            4
## 650
                           13
           990_at
## 651
         991_g_at
                            13
```

To get this in the classic named-list format:

```
z <- as.list(revmap(x)[chroms])
names(z)
## [1] "4" "13"
z[["Y"]]
## NULL</pre>
```

Many of the common methods for accessing *Bimap* objects return things in list format. This can be convenient. But you have to be careful about this if you want to use unlist(). For example the following will return multiple probes for each chromosome:

```
chrs = c("12","6")
mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA)

## $`12`
## [1] "1018_at" "1019_g_at" "101_at" "1021_at"
##
##
##
## $`6`
```

```
## [1] "1026_s_at" "1027_at"
```

But look what happens here if we try to unlist that:

```
unlist(mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA))
## 121 122 123 124 61 62
## "1018_at" "1019_g_at" "101_at" "1021_at" "1026_s_at" "1027_at"
```

Yuck! One trick that will sometimes help is to use Rfunctionunlist2. But be careful here too. Depending on what step comes next, Rfunctionunlist2 may not really help you...

```
unlist2(mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA))
## 12 12 12 12 6 6 6
## "1018_at" "1019_g_at" "101_at" "1021_at" "1026_s_at" "1027_at"
```

Lets ask if the probes in 'pbids' mapped to cytogenetic location "18q11.2"?

To coerce this map to a named vector:

The coercion of the reverse map works too but issues a warning because of the duplicated names for the reasons stated above:

```
cyto2pb <- as.character(revmap(x))
## Warning in .local(x, ...): returned vector has duplicated names</pre>
```

2.0.8 Accessing probes that map to multiple targets

In many probe packages, some probes are known to map to multiple genes. The reasons for this can be biological as happens in the arabidopsis packages, but usually it is due to the fact that the genome builds that chip platforms were based on were less stable than desired. Thus what may have originally been a probe designed to measure one thing can end up measuring many things. Usually you don't want to use probes like this, because if they manufacturer doesn't know what they map to then their usefullness is definitely suspect. For this reason, by default all chip packages will normally hide such probes in the standard mappings. But sometimes you may want access to the answers that the manufacturer says such a probe will map to. In such

cases, you will want to use the toggleProbes method. To use this method, just call it on a standard mapping and copy the result into a new mapping (you cannot alter the original mapping). Then treat the new mapping as you would any other mapping.

```
## How many probes?
dim(hgu95av2ENTREZID)

## [1] 11543    2

## Make a mapping with multiple probes exposed
multi <- toggleProbes(hgu95av2ENTREZID, "all")

## How many probes?
dim(multi)

## [1] 13372    2</pre>
```

If you then decide that you want to make a mapping that has only multiple mappings or you wish to revert one of your maps back to the default state of only showing the single mappings then you can use toggleProbes to switch back and forth.

Finally, there are also a pair of test methods hasMultiProbes and hasSingleProbes that can be used to see what methods a mapping presently has exposed.

```
## Test the multiOnly mapping
hasMultiProbes(multiOnly)

## [1] TRUE
hasSingleProbes(multiOnly)

## [1] FALSE

## Test the singleOnly mapping
hasMultiProbes(singleOnly)

## [1] FALSE
hasSingleProbes(singleOnly)
```

2.0.9 Using SQL to access things directly

While the mapping objects provide a lot of convenience, sometimes there are definite benefits to writing a simple SQL query. But in order to do this, it is necessary to know a few things. The 1st thing you will need to know is some SQL. Fortunately, it is quite easy to learn enough basic SQL to get stuff out of a database. Here are 4 basic SQL things that you may find handy:

First, you need to know about SELECT statements. A simple example would look something like this:

SELECT * FROM genes;

Which would select everything from the genes table.

SELECT gene_id FROM genes;

Will select only the gene_id field from the genes table.

Second you need to know about WHERE clauses:

SELECT gene_id,_id FROM genes WHERE gene_id=1;

Will only get records from the genes table where the gene_id is = 1.

Thirdly, you will want to know about an inner join:

SELECT * FROM genes, chromosomes WHERE genes._id=chromosomes._id;

This is only slightly more complicated to understand. Here we want to get all the records that are in both the 'genes' and 'chromosomes' tables, but we only want ones where the '_id' field is identical. This is known as an inner join because we only want the elements that are in both of these tables with respect to '_id'. There are other kinds of joins that are worth learning about, but most of the time, this is all you will need to do.

Finally, it is worthwhile to learn about the AS keyword which is useful for making long queries easier to read. For the previous example, we could have written it this way to save space:

SELECT * FROM genes AS g,chromosomes AS c WHERE g._id=c._id;

In a simple example like this you might not see a lot of savings from using AS, so lets consider what happens when we want to also specify which fields we want:

SELECT g.gene_id,c.chromosome FROM genes AS g,chromosomes AS c WHERE g._id=c._id;

Now you are most of the way there to being able to query the databases directly. The only other thing you need to know is a little bit about how to access these databases from R. With each package, you will also get a method that will print the schema for its database, you can view this to see what sorts of tables are present etc.

```
org.Hs.eg_dbschema()
```

To access the data in a database, you will need to connect to it. Fortunately, each package will automatically give you a connection object to that database when it loads.

```
org.Hs.eg_dbconn()
```

You can use this connection object like this:

```
query <- "SELECT gene_id FROM genes LIMIT 10;"
result = dbGetQuery(org.Hs.eg_dbconn(), query)
result</pre>
```

Exercise 5

Retrieve the entrez gene ID and chromosome by using a database query. Show how you could do the same thing by using toTable

2.0.10 Combining data from multiple annotation packages at the SQL level

For a more complex example, consider the task of obtaining all gene symbols which are probed on a chip that have at least one GO BP ID annotation with evidence code IMP, IGI, IPI, or IDA. Here is one way to extract this using the environment-based packages:

```
## Obtain SYMBOLS with at least one GO BP
## annotation with evidence IMP, IGI, IPI, or IDA.
system.time({
bpids <- eapply(hgu95av2GO, function(x) {</pre>
    if (length(x) == 1 \&\& is.na(x))
      NA
    else {
         sapply(x, function(z) {
             if (z$Ontology == "BP")
               z$GOID
             else
               NA
             })
})
bpids <- unlist(bpids)</pre>
bpids <- unique(bpids[!is.na(bpids)])</pre>
g2p <- mget(bpids, hgu95av2G02PR0BE)</pre>
wantedp <- lapply(g2p, function(x) {</pre>
    x[names(x) %in% c("IMP", "IGI", "IPI", "IDA")]
})
wantedp <- wantedp[sapply(wantedp, length) > 0]
wantedp <- unique(unlist(wantedp))</pre>
ans <- unlist(mget(wantedp, hgu95av2SYMBOL))</pre>
})
length(ans)
ans[1:10]
```

All of the above code could have been reduced to a single SQL query with the SQLite-based packages. But to put together this query, you would need to look 1st at the schema to know what tables are present:

```
hgu95av2_dbschema()
```

This function will give you an output of all the create table statements that were used to generate the hgu95av2 database. In this case, this is a chip package, so you will also need to see the schema for the organism package

that it depends on. To learn what package it depends on, look at the ORGPKG value:

```
hgu95av20RGPKG
```

Then you can see that schema by looking at its schema method:

```
org.Hs.eg_dbschema()
```

So now we can see that we want to connect the data in the go_bp, and symbol tables from the org.Hs.eg.sqlite database along with the probes data in the hgu95av2.sqlite database. How can we do that?

It turns out that one of the great conveniences of SQLite is that it allows other databases to be 'ATTACHed'. Thus, we can keep our data in many differnt databases, and then 'ATTACH' them to each other in a modular fashion. The databases for a given build have been built together and frozen into a single version specifically to allow this sort of behavoir. To use this feature, the SQLite ATTACH command requires the filename for the database file on your filesystem. Fortunately, R provides a nice system independent way of getting that information. Note that the name of the database is always the same as the name of the package, with the suffix '.sqlite'.:

```
orgDBLoc = system.file("extdata", "org.Hs.eg.sqlite", package="org.Hs.eg.db")
attachSQL = paste("ATTACH '", orgDBLoc, "' AS orgDB;", sep = "")
dbGetQuery(hgu95av2_dbconn(), attachSQL)
```

Finally, you can assemble a cross-db sql query and use the helper function as follows. Note that when we want to refer to tables in the attached database, we have to use the 'orgDB' prefix that we specified in the 'ATTACH' query above.:

```
system.time({
SQL <- "SELECT DISTINCT probe_id,symbol FROM probes, orgDB.gene_info AS gi, orgDB.genes AS g, org
zz <- dbGetQuery(hgu95av2_dbconn(), SQL)
})
##    user    system elapsed
##    0.27    0.06    0.66
#its a good idea to always DETACH your database when you are finished...
dbGetQuery(hgu95av2_dbconn(), "DETACH orgDB"    )</pre>
```

Exercise 6

Retrieve the entrez gene ID, chromosome location information and cytoband information by using a single database query.

Exercise 7

Expand on the example in the text above to combine data from the hgu95av2.db and org.Hs.eg.db with the GO.db package so as to include the GO ID, and term definition in the output.

The version number of R and packages loaded for generating the vignette were:

```
## R version 3.1.3 (2015-03-09)
## Platform: i386-w64-mingw32/i386 (32-bit)
## Running under: Windows Server 2008 R2 x64 (build 7601) Service Pack 1
##
## locale:
```

```
## [1] LC_COLLATE=C
## [2] LC_CTYPE=English_United States.1252
## [3] LC_MONETARY=English_United States.1252
## [4] LC_NUMERIC=C
## [5] LC_TIME=English_United States.1252
##
## attached base packages:
## [1] parallel stats4 stats graphics grDevices utils
## [7] datasets methods base
##
## other attached packages:
## [1] GO.db_3.0.0
                                 hgu95av2.db_3.0.0 AnnotationForge_1.8.2
## [4] org.Hs.eg.db_3.0.0 RSQLite_1.0.0 DBI_0.3.1
## [7] AnnotationDbi_1.28.2 GenomeInfoDb_1.2.4 IRanges_2.0.1
## [10] S4Vectors_0.4.0 Biobase_2.26.0 BiocGenerics_0.12.1
## [13] knitr_1.9
##
## loaded via a namespace (and not attached):
## [1] BiocStyle_1.4.1 evaluate_0.5.5 formatR_1.0 highr_0.4
## [5] stringr_0.6.2 tools_3.1.3
```